

No: 20195-042

## Bipolar forceps Classic, bajonet



### Dimensions and weight

Tip	tip 0.3 mm, very fine
Length	250 mm
Gripping length	6 mm

### Basic features

Shape / Forming	bayonet
-----------------	---------

### Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

### Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

### Additional product information

Max. electrical capacity	0.5 kVp
Packaging unit	1 Piece