

No: 21191-164

## Ball electrode

- The right shape and size for everything and everybody
- Used in all open surgical specialist disciplines
- The entire range available with both 4 mm as well as 2.4 mm shafts



### Dimensions and weight

|                |        |
|----------------|--------|
| Shaft diameter | 2.4 mm |
| Length         | 40 mm  |
| Dimensions     | ø 5 mm |

### Basic features

|                 |               |
|-----------------|---------------|
| Shape / Forming | straight      |
| Compatibility   | VIO, ICC, ACC |
| Flexibility     | rigid         |
| Colour          | black         |

### Cleaning and Reprocessing

|                                 |            |
|---------------------------------|------------|
| Cleaning in washer disinfectant | max 95 °C  |
| Sterilization in autoclave      | max 138 °C |

### Standards

|   |     |
|---|-----|
| Classification acc. to 93/42 EWG        | IIb |
| Protection category acc. to EN 60 601-1 |     |
| Type acc. to EN 60 601-1                |     |

### Additional product information

|                          |         |
|--------------------------|---------|
| Packaging unit           | 5 Piece |
| Max. electrical capacity | 4.3 kVp |