


No: 21191-080

Small surgical set, ø 4 mm



- The right shape and size for everything and everybody
- Used in all open surgical specialist disciplines

Dimensions and weight

Shaft diameter	4 mm
----------------	------

Basic features

Compatibility	VIO, ICC, ACC
Colour	black

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

Additional product information

Packaging unit	1 Piece
Max. electrical capacity	4.3 kVp