

No: 21191-470

Coated Needle Electrode



Dimensions and weight

| | |
|----------------|---------------|
| Dimensions | ø 0,9 x 16 mm |
| Shaft diameter | 2.4 mm |
| Length | 134 mm |

Basic features

| | |
|-----------------|----------|
| Shape / Forming | straight |
| Colour | blue |

Surface

| | |
|---------|---------------------|
| Coating | {PTFE-Beschichtung} |
|---------|---------------------|

Cleaning and Reprocessing

| | |
|--------------------|-----|
| Sterile packaging | E0 |
| Single-use product | yes |

Additional product information

| | |
|--------------------------|----------|
| Packaging unit | 10 Piece |
| Max. electrical capacity | 4 kVp |