

No: 20195-098



Bipolar forceps Classic, straight

Dimensions and weight

Tip	tip 1.1 mm, blunt
Length	190 mm
Gripping length	6 mm

Basic features

Shape / Forming	straight
-----------------	----------

Surface

Branches	corrugated
----------	------------

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 134 °C

Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Max. electrical capacity	0.5 kVp
Packaging unit	1 Piece