

No: 21191-158



## Spatula electrode

- The right shape and size for everything and everybody
- Used in all open surgical specialist disciplines

### Dimensions and weight

Shaft diameter	2.4 mm
Length	180 mm
Dimensions	2,3 x 19 mm

### Basic features

Shape / Forming	straight
Compatibility	VIO, ICC, ACC
Flexibility	rigid
Colour	black

### Surface

Insulation	insulated
------------	-----------

### Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

### Additional product information

Packaging unit	1 Piece
Max. electrical capacity	4.3 kVp