

No: 20196-062

Bipolar connecting cable



Dimensions and weight

Cable length	{4,50 m}
--------------	----------

Basic features

Compatibility	VIO, ICC, ACC
---------------	---------------

Connection	Standard (BI 8 / 4)
------------	---------------------

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
---------------------------------	-----------

Sterilization in autoclave	max 138 °C
----------------------------	------------

Additional product information

Packaging unit	1 Piece
----------------	---------

Max. electrical capacity	1 kVp
--------------------------	-------