

No: 20195-560

## Bipolar forceps PREMIUM, bajonet



### Dimensions and weight

Tip	tip 0,2 mm, pointed
Length	200 mm

### Basic features

Shape / Forming	bayonet
-----------------	---------

### Surface

Branches	roughened
----------	-----------

### Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

### Standards

Classification acc. to 93/42 EWG	IIb
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

### Additional product information

Packaging unit	1 Piece
Max. electrical capacity	0.5 kVp