

No: 20420-021

Cryoprobe tip



Dimensions and weight

| | |
|------------|---------|
| Dimensions | ø 21 mm |
| Length | 18 mm |

Basic features

| | |
|-------------|----------|
| Flexibility | flexible |
|-------------|----------|

Cleaning and Reprocessing

| | |
|---------------------------------|------------|
| Cleaning in washer disinfectant | max 95 °C |
| Sterilization in autoclave | max 134 °C |

Standards

| | |
|---|-----|
| Classification acc. to 93/42 EWG | Ila |
| Protection category acc. to EN 60 601-1 | |
| Type acc. to EN 60 601-1 | |

Additional product information

| | |
|----------------|---------|
| Packaging unit | 1 Piece |
|----------------|---------|