

No: 20420-021

Cryoprobe tip



Dimensions and weight

Dimensions	ø 21 mm
Length	18 mm

Basic features

Flexibility	flexible
-------------	----------

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 134 °C

Standards

Classification acc. to 93/42 EWG	Ila
Protection category acc. to EN 60 601-1	
Type acc. to EN 60 601-1	

Additional product information

Packaging unit	1 Piece
----------------	---------