

No: 20194-077

Patient plate cable



Dimensions and weight

Cable length	{4 m}
--------------	-------

Basic features

Compatibility	VIO, ICC, ACC, T-Series
---------------	-------------------------

Connection	Standard
------------	----------

Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
---------------------------------	-----------

Sterilization in autoclave	max 138 °C
----------------------------	------------

Additional product information

Packaging unit	1 Piece
----------------	---------

Max. electrical capacity	0.5 kVp
--------------------------	---------