

No: 21191-012

## Needle electrode

- The right shape and size for everything and everybody
- Used in all open surgical specialist disciplines



### Dimensions and weight

Shaft diameter	4 mm
Length	40 mm
Dimensions	ø 0,8 x 22 mm

### Basic features

Shape / Forming	straight
Compatibility	VIO, ICC, ACC
Flexibility	rigid
Colour	black

### Cleaning and Reprocessing

Cleaning in washer disinfectant	max 95 °C
Sterilization in autoclave	max 138 °C

### Additional product information

Packaging unit	5 Piece
Max. electrical capacity	4.3 kVp